

Name \_\_\_\_\_ Contact Name \_\_\_\_\_  
 Bill To \_\_\_\_\_ City \_\_\_\_\_ State \_\_\_\_\_ Zip \_\_\_\_\_  
 Ship To (if different) \_\_\_\_\_ City \_\_\_\_\_ State \_\_\_\_\_ Zip \_\_\_\_\_  
 Phone \_\_\_\_\_ Fax \_\_\_\_\_ Email \_\_\_\_\_  
 Visa, MC, Discover \_\_\_\_\_ Expiration Date \_\_\_\_\_

**For multiple case orders of developer and fixer - indicate the TOTAL number of each product. Thanks!**

<b>AUTOMATIC PROCESSING (ROLLER &amp; PERI-PRO*)</b>	<b>CASE PACK</b>	<b>ITEM #</b>	<b>QTY.</b>	<b>PRICE</b>	<b>TOTAL</b>
<b>Tech Chem™ Dev / Fixer Premix 4 Gallons per case (Roller)</b>	Dev + Fix	TCP4128		\$ 36.00	
<b>Peri-Kem® Dev / Fixer Premix 6 Quarts (Peri-Pro machines*)</b>	Dev + Fix	PKP6032		\$ 26.00	
<b>Peri-Kem® Dev / Fixer Premix 12 Quarts (Peri-Pro*)</b>	Dev + Fix	PKP1232		\$ 48.00	
<b>Dental Processor Cleaner™ 6 cleanings per case</b>	1 Case	DPC0664		\$ 46.00	
<b>Dental Processor Cleaner Bulk up to 8 cleanings</b>	1 lb 3 oz Jar	DPC1001		\$ 32.00	
<b>Dental Processor Cleaner Super Saver up to 42 cleanings</b>	6 lb Jar	DPC1005		\$ 120.00	

<b>MANUAL PROCESSING (Dip &amp; Dunk Tanks)</b>	<b>CASE PACK</b>	<b>ITEM #</b>	<b>QTY.</b>	<b>PRICE</b>	<b>TOTAL</b>
<b>XLR™ Developer / Fixer Concentrate 6 Quarts per case</b>	Dev + Fix	XLR6032		\$ 38.00	
<b>XLR™ Developer / Fixer Concentrate 12 Quarts per case</b>	Dev + Fix	XLR1232		\$ 68.00	
<b>Dip Tank Cleaner™ 4 cleanings</b>	1 lb 3 oz Jar	DTC1216		\$ 32.00	
<b>Dip Tank Cleaner Super Saver 21 cleanings</b>	6 lb Jar	DTC1006		\$ 120.00	

<b>MISCELLANEOUS</b>	<b>UNIT</b>	<b>ITEM #</b>	<b>QTY.</b>	<b>PRICE</b>	<b>TOTAL</b>
<b>Fixer Check™ for X-Ray Use up to a 3-year supply</b>	1 oz Bottle	FOC01001		\$ 10.00	
<b>Stain Away™ Clothing Stain Remover</b>	12 oz Bottle	CSR01012		\$ 14.00	
<b>VacAttack™ Enzyme Cleaner Conc. w/ Pump Dispenser</b>	1 Gallon	VAC1128		\$ 30.00	
<b>VacAttack™ Enzyme Cleaner Conc. w/ 1 Pump Dispenser</b>	4 Gallons	VAC4128		\$100.00	

\*Not a registered trade mark of Tech Chem, Inc.

Subtotal

Missouri Residents add 6.475% Tax

**Total**

Prices subject to change without notice. **ACTUAL shipping will be added to your order! No handling fees - No minimums!!**

**Order securely on line [www.techcheminc.com](http://www.techcheminc.com), Mail, Fax, Email, or Phone your completed order to:**

Tech Chem, Inc.  
 321 Route Y  
 Arcola MO 65603-7703

**Toll-Free Phone 800-852-5641 • sales@techcheminc.com • Toll-Free Fax 888-424-3294** Rev 12/10

# TECH CHEM INC

*Where your image matters.*

**Tech Chem's Developer & Fixer Premix™** is specially formulated to provide sharp, crisp x-rays in ALL Automatic Film Processors that have ROLLER TYPE transport assemblies. These include AT2000®, Dent-X 810® & 9000®, Velopex Intra®, Velopex Xtra®, and Velopex Xtender® processors to name a few.\* This is packaged in cases of 4 gallons of ready-to-use solution. You may order any combination of developer and/or fixer in full cases to suit your needs. **Without paying extra!**

**4 Gallons per case \$ 36.00**

**Peri-Kem®** is manufactured to produce superior x-rays in ALL Automatic Film Processors that have ROLLERLESS Type transport assemblies, like the Peri-Pro® machine.\*

This is packaged in cases of 6 32-ounce bottles OR 12 32-ounce bottles of ready-to-use solution. You may order any combination of developer and / or fixer in full cases to suit your needs. **Without paying extra!**

**6 32-oz Bottles for \$ 26.00 - 12 32-oz Bottles for \$ 48.00**

**Dental Processor Cleaner™** is designed to clean ALL Automatic Film Processors - roller & rollerless! Simply activate with HOT Tap Water, wait 20-30 minutes, drain and rinse thoroughly, replace transports & you're ready to go!

**No hazardous chemicals! Biodegradable! Safe to YOU!**

**This is packaged 3 ways:** Pre-measured in six convenient bottles - **6 cleanings for \$46.00;**

Compactly in a 1# jar - **up to 8 cleanings (scoop included) for \$32.00;**

Economically in a 6# canister - **up to 42 cleanings (scoop included) for \$120.00**

**Fixer Check for X-Ray Use™** is for use in determining Silver Saturation in fixer for ALL types of fixer - Automatic or Manual! Using the dropper provided, put 2 drops of Fixer Check directly in the fixer. If the spot turns milky white your fixer has too much silver in it. If the spot remains clear, YOUR FIXER IS FINE & STILL USABLE!

**Fixer outlasts developer 2-6 times!** Stop throwing away - or worse paying for hauling away - good fixer! One bottle of Fixer Check for X-Ray Use will last **UP TO 3 YEARS!** You may purchase developer and fixer from us the way you NEED IT! Without paying extra!

**1 oz Bottle for \$ 10.00 - UP TO A 3 YEAR SUPPLY!**

**XLR™ Developer and Fixer Concentrate** is formulated to provide sharp, crisp x-rays in all Manual Tanks (Dip & Dunk). Each bottle of XLR concentrate makes 1 Gallon of working solution!

This is packaged in cases of 6 32-ounce bottles OR 12 32-ounce bottles of concentrate solution. You may order any combination of developer and / or fixer in full cases to suit your needs. Without paying extra!

**6 32-oz Bottles for \$ 38.00 - 12 32-oz Bottles for \$ 68.00**

**Dip Tank Cleaner™** is activated with HOT Tap Water. Drain chemicals from the tanks, add Dip Tank Cleaner, Hot Tap Water, wait 20-30 minutes, drain, rinse & you're done! No more scrubbing!

This is packaged 2 ways: Compactly in a 1# jar - **4 cleanings (scoop included) for \$32.00;** Economically in a 6# canister - **21 cleanings (scoop included) for \$120.00**

**VacAttack™ Enzymatic Cleaner Concentrate** is for use in ultrasonic units, evacuation systems, collection/holding tanks, and as an instrument presoak. One gallon of VacAttack makes up to 42 gallons of ready-to-use solution. Why pay more? Completely removes any organic substance, including dried-on blood and tissue from rubber, glass, plastics, all metals (including aluminum), and fabric.

Enzymatic cleaner tips: DO NOT use any enzymatic cleaner (regardless of the manufacturer) at temperatures above 170°F (76°C) as the enzyme will be destroyed. DO NOT combine with bleach, detergents, or any other chemicals. Sterilization or disinfecting require a separate, subsequent step.

**1 Gallon for \$ 30.00 - 4 Gallons for \$ 100.00**

**Stain Away™** removes most stains caused by developer and fixer from clothing and other fabrics! Safely remove stains caused by Developer and Fixer from clothing and other fabrics. This product has been known to even remove stains that were several years old. (Removing old developer stains from athletic shoes is a favorite!)

**12 oz Bottle for \$ 14.00**

***Actual shipping charges are added to each order - no handling fees!!***

\*Not registered trademarks of Tech Chem, Inc.